THIS PAGE IS INSERTED BY OIPE SCANNING AND IS NOT PART OF THE OFFICIAL RECORD

Best Available Images

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

BLACK BORDERS

TEXT CUT OFF AT TOP, BOTTOM OR SIDES

FADED TEXT

BLURRY OR ILLEGIBLE TEXT

SKEWED/SLANTED IMAGES

COLORED PHOTOS HAVE BEEN RENDERED INTO BLACK AND WHITE

VERY DARK BLACK AND WHITE PHOTOS

UNDECIPHERABLE GRAY SCALE DOCUMENTS

IMAGES ARE THE BEST AVAILABLE COPY. AS RESCANNING WILL NOT CORRECT IMAGES, PLEASE DO NOT REPORT THE IMAGES TO THE PROBLEM IMAGE BOX.






Figure 1



Approximate flow for straight tubes







FIGURE 8

Approximate flow for wiggle (variable width tubes)








FIGURE 9A



pressure line pressure is.

7 individual high pressure




Bake ~2 hrs 65°C


Seal to microscope slide

Cut out individual devices

and peel from mold

Another variation:






Seal two together back to back for double sized chamber.

Another embodiment




Same elastomer molded tubes, or etched glass tubes, or laser machined tubes in plastic





Conventionally machined.
wells, or injection molded
wells, or epoxy molded
wells, or elatomer molded
wells.



Thin 100-300um thick PDMS membrane spin on and peel off


FIGURE 12

Well or mixing chamber Injection port (4) And channel

each well is 5mm x 10mm. The spacing between the This is the mask for the The dashed circle is 3.5 well, or chamber mold. inches in diameter and wells is 1mm. The dark areas on the mask areas will be cavities in the After casting the elastomer substrate to form the mold. 500µm) on a silicon wafer and peeling it off, the dark structures (typ: 300will result in raised silicone elastomer. This is the mask for three different kinds of tubes.

The left is the "wiggle tube" to drive flow transverse to the peristalsis direction. The middle is the standard straight tube, and the right are straight tubes with narrower regions near the edges of the wells. (The next slide shows the two masks overlayed)

The dark areas on the mask will result in raised structures (typ: 200-400µm) on a silicon wafer substrate to form the mold. After casting the elastomer and peeling it off, the dark areas will be cavities in the silicone elastomer.




This shows the two masks overlayed.


A key part of the design is the scalability. Note that only seven macroscopic actuators (Lee company solenoid valves – see later slide). Are required to actuate 8 wells because each tube crosses over eight wells. By making it longer, scalability would improve. (7 valves actuate 16, or 32 wells)


Current implementations are limited by the clean room processing tubes, and easy availability of 1x3 inch microscope slides substrates.

Minimizing the number of actuators is important to reduce overall system size and to keep to simplify the assembly of the macroscopic actuators.



Exploded view of manifold, which mainly serves to reduce the pitch of the solenoid valves down to the pitch of the acutation ports.









FIGURE 20



Note, the light areas look light because the tubes are deflected and the dye is pushed away, allowing more white background to show through. 8psi, 100Hz.






FIGURE 22

Mixing Results



Slow, 10Hz actuation rate, 8psi Corn syrup



more separated flow, symmetry broken, no vortices don't mix Upper/lower half but sometimes



Corn syrup

Fast, 100Hz actuation rate, 8psi






∆t=8sec



∆t=16sec



∆t=32sec

Mixing trends

